DEVILBISS

515 Series
Serie de 515
Série 515
515 Serie
Serie da 515

 **SUNRISE
MEDICAL**

DeVilbiss[®] Oxygen Concentrator

Instruction Guide

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TABLE OF CONTENTS

Specifications.....	3
Important Safeguards	9
Introduction.....	9
Why Your Physician Prescribed Supplemental Oxygen	10
How Your Concentrator Works	10
Important Parts of Your Concentrator.....	11
Setting Up Your Concentrator.....	13
Before Operating Your Concentrator.....	13
Operating Your Concentrator.....	15
Reserve Oxygen System	17
Caring for Your Concentrator.....	17
Minor Troubleshooting	19
Minor OSD® Troubleshooting	19

CAUTION—Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

INDICATIONS FOR USE— The DeVilbiss Oxygen Concentrator is intended for use as an oxygen concentrator to provide supplemental low flow oxygen therapy in the home, nursing homes, patient care facilities, etc.

⚠ WARNING—Under certain circumstances, oxygen therapy can be hazardous. Seeking medical advice before using an oxygen concentrator is advisable.

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SPECIFICATIONS

	DEVILBISS 4-LITER SERIES	DEVILBISS 5-LITER SERIES	
Catalog Number	515UK	515DS, 515DZ	515KS, 515KZ, 515NS
Delivery Rate (Lower delivery rates available for low flow applications)	1 to 4 LPM	1 to 5 LPM	
Maximum Recommended Flow	4 LPM	5 LPM	
Outlet Pressure	8.5 psig (58.6 kPa)		
Electrical Rating	230 V, 50 Hz, 1.4 Amp	115 V, 60 Hz, 4.2 Amp	220-230 V~, 50 Hz, 1.7 Amp 230 V~, 60 Hz, 1.9 Amp
Operating Voltage Range	195-253 V~, 50 Hz	97-127 V~, 60 Hz	187-255 V~, 50 Hz 218-253 V~, 60 Hz
Oxygen Percentage	1-4 LPM = 93%±3%	1-5 LPM=93%±3%	
Operating Altitude (tested at 21°C only) 0-1500 M (0-4921 ft) 1500-3000 M (4921-9842 ft) 3000-4000 M (9842-13123 ft)	Across the voltage range: No degradation of performance Not recommended/Not tested Not recommended/ Not tested	Across the voltage range: No degradation of performance Tested at nominal voltage only: No degradation in performance Not recommended/Not tested	Across the voltage range: No degradation in performance Tested at 230V/50Hz only: ~4% oxygen performance drop at 5 LPM Not recommended/Not tested
Operating Environment Range* 10°C to 35°C, humidity range of 30% to 75%	No degradation in performance across the operating voltage range (tested at sea level).	No degradation in performance across the operating voltage range (tested at sea level).	~4% oxygen performance drop at 5 LPM and 35°C condition. No degradation in performance across the operating voltage range at other liter flows (tested at 670 M).
Power Consumption	295 Watts Average	400 Watts Average	230v / 50 Hz - 365 Watts Average 230v / 60 Hz - 435 Watts Average
Weight	53 lbs. (24.5 Kilograms)	52 lbs. (23.5 Kilograms)	53 lbs. (24.5 Kilograms)
Sound Level (ISO 8359:1996 from front)	50.5 dbA overall average	52.5 dbA overall average**	50.5 dbA overall average (50Hz)
Dimensions	27.75"H x 16"W x 14"D (70.5 x 40.6 x 35.6 cm)		
Pressure Relief Valve	44 psig±3psig (303 kPa±21 kPa)		
Operating System	Time Cycle / Pressure Swing		
For units equipped with an OSD, the visible "low oxygen" indicator will activate at the following level	Units with serial numbers less than H20000: 85%±2% (The audible alarm will alert at approximately 75%) Units with serial numbers of H20000 and greater: 83.5%±1.5% (The audible alarm will alert at approximately 75%. At less than 60%, the red "service required" light will activate.)		
Storage Conditions	-40°C to 70°C, humidity range of 10% to 100%, including condensation		
Equipment Class and Type	<input type="checkbox"/> Class II Equipment Double Insulated <input checked="" type="checkbox"/> Type B Applied Part		
Approval Body and Safety Standard	TUV IEC 601-1+A1+A2 EN 60601+A1+A2 ISO 8359:1996	CSA CAN/CSA-C22.2 No. 601.1-M90	TUV IEC 601-1+A1+A2 EN 60601+A1+A2 ISO 8359:1996 (515KZ excluded)
EMC Compliance To:	IEC 601-1-2		

*Note: The OSD performance at 10°C to 35°C, 95% R.H. through voltage range on the 515DS verified at 670m.

**Note: This corresponds with a sound level of 50 dbA overall average as measured according to the guidelines of the now defunct ANSI Z79.13-1981 standard.

Specifications subject to change without notice.

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DECLARATION OF CONFORMITY**Manufacturer:**

Sunrise Medical
Respiratory Products Division

Address:

100 DeVilbiss Drive
Somerset, PA 15501 USA

Product Designation: Oxygen Concentrators

Catalog Number:

515UK, 515KS, 515NS

We herewith declare that the above-mentioned product complies with the requirements of EC directive 93/42/EEC and the following (see Specifications Matrix for specific compliance by model):

Quality System Standards Applied:

ISO9001/EN46001

Notified body RWTÜV**MDD Annex II Applied****Safety Standards Applied:**

IEC 601-1+A1+A2
EN 60601+A1+A2
ISO8359: 1996

EMC Compliance To:

IEC 601-1-2

Place: Somerset, PA

Signature:

Name: Michael Chellson

Title: Director of quality and regulatory affairs

European Contact:

Sunrise Medical Ltd.
Sunrise Business Park
High Street
Wollaston, West Midlands
DY8 4PS
ENGLAND
44-138-444-6688

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Nome: Michael Chellson

Titolo:

Director of quality and regulatory affairs

Contatto in Europa:

Sunrise Medical Ltd
Sunrise Business Park
High Street
Wollaston, West Midlands
DY8 4PS
ENGLAND
44-138-444-6688

0044

IMPORTANT SAFEGUARDS

⚠ Read this entire guide before using your DeVilbiss concentrator. Important information is highlighted by these terms throughout this guide:

DANGER—Urgent safety information for hazards that will cause serious injury or death.

WARNING—Important safety information for hazards that might cause serious injury.

CAUTION—Information for preventing damage to the product.

NOTE—Information to which you should pay special attention.

Important safeguards are indicated throughout this guide; pay special attention to all safety information.

INTRODUCTION

This instruction guide will acquaint you with the DeVilbiss oxygen concentrator. Make sure that you read and understand this guide before operating your unit. Important safeguards are indicated throughout this guide; pay special attention to all safety information. Contact your Sunrise Medical equipment provider should you have any questions.

⚠ WARNING—For your safety, the oxygen concentrator must be used according to the prescription determined by your physician.

⚠ DANGER—Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy. Do not use within 5 feet (1.6m) of hot, sparking objects or naked sources of flame.

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Why Your Physician Prescribed Supplemental Oxygen

Today, many people suffer from heart, lung, and other respiratory diseases. Many of these people can benefit from supplemental oxygen therapy. Your body requires a steady supply of oxygen to function properly. Your physician prescribed supplemental oxygen for you because you are not getting enough oxygen from room air alone. Supplemental oxygen will increase the amount of oxygen that your body receives.

Supplemental oxygen is not addictive. Your physician prescribed a specific oxygen flow to improve symptoms such as headaches, drowsiness, confusion, fatigue, or increased irritability. If these symptoms persist after you begin your supplemental oxygen program, consult your physician.

How Your DeVilbiss Oxygen Concentrator Works

Oxygen concentrators are the most reliable, efficient, and convenient source of supplemental oxygen available today. The oxygen concentrator is electrically operated. The unit separates oxygen from room air which allows high-purity supplemental oxygen to be delivered to you through the oxygen outlet. Although the concentrator filters the oxygen in a room, it will not affect the normal amount of oxygen in your room.

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IMPORTANT PARTS OF YOUR CONCENTRATOR

Please take time to familiarize yourself with your DeVilbiss oxygen concentrator before operating.

Front View

1. Operating instructions
2. Green Power light – illuminates when your concentrator is operating.
3. Power Switch
I = ON
O = OFF
4. Flow meter knob
5. Flow meter
6. Circuit breaker – (selected units) resets the unit after electrical overload shutdown.
7. Oxygen outlet – oxygen is dispersed through this port.
8. Normal Oxygen (green) light – ONLY on units with the OSD® (see page 15).
9. Low Oxygen (yellow) light – ONLY on units with the OSD® (see page 15).
10. Red Service Required light – when illuminated contact your Sunrise Medical provider.

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
Figure 2.
Figura 2.
Abb 2.

Back View

11. Important Safeguards
12. Handgrip
13. Exhaust
14. Power cord and/or IEC power connector.
15. Line cord strap
16. Air filter – prevents dirt, dust, and lint from entering your unit.

Accessories

The accessories below are approved for use with DeVilbiss Oxygen Concentrators:

- Bubble Humidifier (200 ml Capacity)
- 50/caseHUM16
- Cannula (disposable)
- 50/caseCAN00
- Cannula w/ 2.1 m (7 ft)
Oxygen Tubing (disposable)
- 50/caseCAN70
- Oxygen Outlet Connector
- Plastic - 1/packCN100
- Metal - 1/packCN200
- Oxygen Tubing (disposable)
- 2.1 m (7 ft) - 50/caseOST07
- 4.2 m (14 ft) - 50/caseOST15
- 7.5 m (25 ft) - 25/caseOST25
- 15 m (50 ft) - 15/caseOST50

NOTE—The use of certain humidifiers and accessories not specified above may impair the oxygen concentrator's performance.

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SETTING UP YOUR OXYGEN CONCENTRATOR

1. Position your unit near an electrical outlet in the room where you spend most of your time.

⚠ DANGER—Keep the oxygen concentrator at least 5 feet (1.6 m) from hot, sparking objects or naked sources of flame.

NOTE—Do not connect to an electrical outlet controlled by a wall switch. No other appliances should be plugged into the wall outlet.

2. Position your unit at least 6 inches (16 cm) from walls, draperies, or any other objects that might prevent the proper flow of air in and out of your oxygen concentrator. The oxygen concentrator should be located so as to avoid pollutants or fumes.

BEFORE OPERATING YOUR OXYGEN CONCENTRATOR

1. Before operating your unit, always check to be sure the air filter (located on the back of your unit) is clean. Proper cleaning of this filter is discussed in the Caring For Your Concentrator section on page 17.
2. Attach the appropriate oxygen accessories to the oxygen outlet.

Oxygen Tubing Connection:

- a) Thread the oxygen outlet connector onto the oxygen outlet.
- b) Attach the oxygen tubing directly to the connector. **Figure 3.**

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Figure 4.
Figura 4.
Abb 4.

Figure 5.
Figura 5.
Abb 5.

Oxygen Tubing Connection With Humidification:

If your physician has prescribed an oxygen humidifier as part of your therapy, follow these steps:

- Fill the humidifier bottle with distilled water. Do not overfill. (If using a prefill, go to letter b.)
 - Thread the wing nut located on the top of the humidifier bottle to the oxygen outlet so that it is suspended **Figure 4**. Make sure it is securely tightened.
 - Attach the oxygen tubing directly to the humidifier bottle outlet fitting **Figure 5**.
- Your physician has prescribed either a nasal cannula, catheter, or face mask. In most cases, they are already attached to the oxygen tubing. If not, follow the manufacturer's instructions for attachment.
 - Remove the power cord completely from the line cord strap. Make sure the power switch is in the "Off" position, and insert the plug into the wall outlet. The unit is double insulated to guard against electric shock.

NOTE—(only 115 volt units) The plug on the DeVilbiss oxygen concentrator has one blade wider than the other. To reduce the risk of electric shock, this plug is intended to fit in a wall outlet only one way. Do not attempt to defeat this safety feature.

⚠ WARNING—Improper use of the power cord and plugs can cause a burn, fire, or other electric shock hazards. Do not use the unit if the power cord is damaged.

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OPERATING YOUR DEVILBISS OXYGEN CONCENTRATOR

⚠ DANGER—Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy. Keep the oxygen concentrator at least 5 feet (1.6 m) from hot, sparking objects or naked sources of flame.

1. Press the power switch to the "On" position. When the unit is turned on, the Power light will illuminate, and the patient alert system – Service Required light and audible signal – will briefly alarm.

DEVILBISS OXYGEN CONCENTRATORS WITH OSD ONLY

The OSD (Oxygen Sensing Device) is a device within your concentrator that monitors the oxygen produced by your unit. When the unit is turned "On," all four lights (Power, Service Required, Low Oxygen, and Normal Oxygen) on the front panel will illuminate briefly. After a few seconds, only the Power and Normal Oxygen lights will remain lit.

The OSD lights on the front panel are defined as follows:

- Green Normal Oxygen light—acceptable oxygen level.
- Yellow Low Oxygen light—below an acceptable oxygen level.

If the oxygen purity falls below the acceptable level, the green Normal Oxygen light will shut off and the yellow Low Oxygen light will illuminate. Switch to your reserve oxygen system. Refer to the Minor OSD Troubleshooting section in this guide on page 19 and contact your Sunrise Medical provider.

As an added safety feature should the oxygen purity continue to drop, an intermittent audible signal will sound. Contact your Sunrise Medical provider immediately. Do not attempt any other maintenance.

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Figure 6.
Figura 6.
Abb 6.

NOTE– If the Service Required light illuminates and the audible signal alarms but the unit is not operating, there is no power to the unit. Refer to the Minor Troubleshooting chart on page 19 and contact your Sunrise Medical provider if necessary.

2. Check the flow meter to make sure that the flow meter ball is centered on the line next to the prescribed number of your flow rate.

CAUTION–It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen– consult your physician.

NOTE–Your Sunrise Medical provider may have preset the flow meter so that it can not be adjusted.

NOTE–If the flow meter knob is turned clockwise, the flow decreases (and eventually will shut off the oxygen flow). If the knob is turned counterclockwise, the flow increases.

NOTE–On later models of the 515 concentrators equipped with an OSD, the low-flow alarm may be activated if the flow meter ball is set below 1 lpm. The unit will continue to run; however, the Service Required light will come on accompanied by an audible alarm. Adjust the flow meter to your prescribed flow.

3. Your DeVilbiss concentrator is now ready for use, properly position the cannula or mask **Figure 6**. Allow 20 minutes for oxygen concentrator to reach stated performance.

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RESERVE OXYGEN SYSTEM

As a precaution, your Sunrise Medical provider should supply you with a reserve oxygen system. If your unit loses electrical power or fails to operate correctly, the Patient Alert System will sound to signal you to switch to your reserve oxygen system (if provided) and contact your Sunrise Medical provider.

CARING FOR YOUR DEVLBISS OXYGEN CONCENTRATOR

NOTE—Use no lubricants, oils or grease.

⚠ WARNING—Before attempting any cleaning procedures, turn the unit “Off.”

Oxygen Humidifier (reusable bottles)

If your physician has prescribed humidification, clean your humidifier bottle daily. Follow the instructions supplied by the manufacturer. If no cleaning instructions were supplied, follow these steps:

1. Wash the humidifier in a solution of hot water and dish-washing detergent **Figure 7**.
2. Soak the humidifier in a solution of one part white vinegar to three parts hot water for 30 minutes **Figure 8**. This solution acts as a germicidal agent.
3. Rinse thoroughly with hot tap water and refill with distilled water **Figure 9** for use. Do not overfill.

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Figure 9.
Figura 9.
Abb 9.

Cannula, Mask, and Tubing

Clean and replace the cannula, mask, and tubing as instructed by your Sunrise Medical provider.

Air Filter and Oxygen Outlet Connector

The air filter and connector should be cleaned at least once a week.

To clean, follow these steps:

1. Remove the air filter, located on the back of the unit. Remove the oxygen outlet connector (if used).
2. Wash in a solution of warm water and dishwashing detergent **Figure 10**.
3. Rinse thoroughly with warm tap water and towel dry. The filter should be completely dry before reinstalling.

CAUTION—To prevent product damage, do not attempt to operate the unit without the air filter or while the filter is still damp.

Exterior Cabinet

⚠ WARNING—To avoid electric shock, do not remove the concentrator cabinet. The cabinet should only be removed by a qualified Sunrise Medical technician. Do not apply liquid directly to the cabinet or utilize any petroleum-based solvents or cleaning agents.

Clean the concentrator exterior cabinet by using a damp cloth or sponge with a mild household cleaner and wipe it dry.

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
Figure 10.
Figura 10.
Abb 10.

MINOR TROUBLESHOOTING

The following troubleshooting chart will help you analyze and correct minor oxygen concentrator malfunctions. If the suggested procedures do not help, switch to your reserve oxygen system and call your Sunrise Medical homecare provider. Do not attempt any other maintenance.

⚠ WARNING—To avoid electric shock hazard, do not remove the cabinet. The cabinet should only be removed by a qualified Sunrise Medical homecare technician.

MINOR TROUBLESHOOTING CHART

Symptom	Possible Cause	Remedy
A. Unit does not operate. Power light is off when the power switch is "On." Audible alarm is pulsing and Service Required light is flashing.	<ol style="list-style-type: none"> 1. Power cord not properly inserted into wall outlet. 2. No power at wall outlet. 3. Oxygen concentrator circuit breaker activated (selected units). 	<ol style="list-style-type: none"> 1. Check power cord connection at the wall outlet. On 230 volt units, also check the back of the unit. 2. Check your home circuit breaker and reset if necessary. Use a different wall outlet if the situation occurs again. 3. Press the concentrator circuit breaker reset button (if equipped) located below the power switch. Use a different wall outlet if the situation occurs again. <p>If the above remedies do not work, contact your Sunrise Medical provider.</p>
B. Unit operates, the Power light is on when the power switch is "On." Red Service Required light is illuminated. Audible alarm may be sounding.	<ol style="list-style-type: none"> 1. Air filter is blocked. 2. Exhaust is blocked. 3. Blocked or defective cannula, catheter, face mask, or oxygen tubing. 4. Blocked or defective humidifier bottle. 5. Flow meter set too low. 	<ol style="list-style-type: none"> 1. Check the air filter. If the filter is dirty, wash it following the cleaning instructions on page 18. 2. Check the exhaust area; make sure there is nothing restricting the unit exhaust. 3. Detach cannula, catheter, or face mask. If proper flow is restored, clean or replace if necessary. Disconnect the oxygen tubing at the oxygen outlet. If proper flow is restored, check oxygen tubing for obstructions or kinks. Replace if necessary. 4. Detach the humidifier from the oxygen outlet. If proper flow is obtained, clean or replace humidifier. 5. Set flow meter to prescribed flow rate. <p>If the above remedies do not work, contact your Sunrise Medical provider.</p>
C. If any other problems occur with your oxygen concentrator.		<ol style="list-style-type: none"> 1. Turn your unit "Off." Switch to your reserve oxygen system, and contact your Sunrise Medical provider immediately.

MINOR TROUBLESHOOTING CHART - ONLY OSD Concentrators

Visual Indicator Lights	Possible Cause	Remedy
A. Both the green Normal Oxygen and the yellow Low Oxygen lights are either on or off.	<ol style="list-style-type: none"> 1. OSD malfunction. 	<ol style="list-style-type: none"> 1. Contact your Sunrise Medical provider.
B. Yellow Low Oxygen light is on or the yellow Low Oxygen light is on and the intermittent audible signal is sounding.	<ol style="list-style-type: none"> 1. Flow meter is not properly set. 2. Air filter is blocked. 3. Exhaust is blocked. 	<ol style="list-style-type: none"> 1. Ensure the flow meter is properly set to the prescribed number. 2. Check the air filter. If the filter is dirty, wash it following the cleaning instructions on page 18. 3. Check the exhaust area; make sure there is nothing restricting the unit exhaust.
C. Red Service Required light is on and an intermittent audible signal is sounding.	<ol style="list-style-type: none"> 1. Flow meter is not properly set. 2. Air filter is blocked. 3. Exhaust is blocked. 	<ol style="list-style-type: none"> 1. Ensure the flow meter is properly set to the prescribed number. 2. Check the air filter. If the filter is dirty, wash it following the cleaning instructions on page 18. 3. Check the exhaust area; make sure there is nothing restricting the unit exhaust. <p>If the above remedies do not work, contact your Sunrise Medical provider.</p>

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Sunrise Medical Canada, Inc.
Sunrise Medical Pty. Limited